

OVERALL GOALS

- ▶ Creation of a plasmid containing *FAR1-22/FAR1-WT* driven by *GALpr*.
- ▶ Testing under a galactose and glucose gradient to determine level of activation required to stop cell cycle.
- ▶ Flow cytometry to measure DNA content after cell arrest.
- ▶ Creation of a genomic strain containing *FAR1-22/FAR1-WT* driven by *ZEV₂₆₈pr*.
- ▶ Testing under a β -estradiol gradient mimicking the gal/glu gradient.
- ▶ Determination of appropriate conditions for both β -estradiol and blue light induction.
- ▶ Blue light controls.

STEPS FOR PLASMID CONSTRUCTION

- ☒ Miniprep and digest pMM86
- ☒ Miniprep pGAL-FAR1-22 (pTCN113)
- ☒ PCR amplify FAR1-22 fragment
- ☒ Colony PCR amplify FAR1-WT fragment
- ☐ Digest both fragments
- ☐ Gel extraction of fragments
- ☐ Ligation
- ☐ Transformation into bacteria
- ☐ Miniprep and sending to sequencing
- ☐ Transformation into yeast

STEPS FOR GENOMIC STRAIN CONSTRUCTION

- ☒ Miniprep pMM299 (ZEVpr)
- ☒ Miniprep pCORE
- ☒ Miniprep pGAL-FAR1-22 (pTCN113)
- ☒ PCR amplify CORE fragment
- ☒ PCR amplify FAR1-22 fragment
- ☒ Colony PCR amplify FAR1-WT fragment
- ☒ Colony PCR amplify ZEVpr fragment
- ☐ Transform ZEV yeast with CORE
- ☐ Transform ZEV+CORE with ZEVpr + FAR fragment
- ☐ Send to sequencing.
- ☐ Transform ZEV + FAR with the blue light plasmids

OTHER STEPS IN THE NEAR FUTURE:

- ☐ Order sequencing primers.
- ☐ Build blue light rig.
- ☐ Transform ZEV strain w/ Blue light plasmids.
- ☐ Set up Galactose+Glucose stock.
- ☐ Plan for FACS experiment.